

APPLICATION FOR OSHPD SPECIAL SEISMIC	OFFICE	USE ONLY
CERTIFICATION PREAPPROVAL (OSP)	APPLICATION #:	OSP – 0493 – 10
OSHPD Special Seismic Certification Preapproval (OSP)		
Type: 🛛 New 🗌 Renewal		
Manufacturer Information		
Manufacturer: Siemens Healthcare GmbH		
Manufacturer's Technical Representative: Dr. Damian Kopyto		
Mailing Address: Siemensstr. 3, D-91301 Forchheim, Germany		
Telephone: +499191 188778 Email: damiar	n.kopyto@siemens.com	
Product Information		
Product Name: Multitom Rax System		
Product Type: Robotic X-ray medical imaging system		
Product Model Number:       See Attachment         (List all unique product identification numbers and/or part numbers)         General Description:       Multiple component system for producing X-Radiagnostic results         Mounting Description:       Multiple – Mix of rigid floor mounting, combined		
See attachment.		
Applicant Information		
Applicant Company Name: <u>W.E. Gundy &amp; Associates, Inc.</u>		
Contact Person: Travis Soppe, SE		
Mailing Address: _250 Bobwhite Ct, Suite 100, Boise, ID 83706		
Telephone: (208) 342-5898 Ext. 115 Email: tsoppe	@wegai.com	
I hereby agree to reimburse the Office of Statewide Health I accordance with the California Administrative Code, 2016.	Planning and Develo	pment review fees in
Signature of Applicant:	Date:	12-08-2016
Title: Vice President Company Name: W.E. C	Sundy & Associates, Inc.	
"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"	All All Ann	OSHPD
STATE OF CALIFORNIA – HEALTH AND HUMAN SERVICES AGENCY OSH-FD-759 (REV 12/16/15)	JMAAAAA	Page 1 of 3



California Licensed Structural Engineer Responsible for the Engineering and Test Report(s)
Company Name: W.E. Gundy & Associates, Inc.
Name: Travis Soppe, SE California License Number: S6115
Mailing Address: 205 Bobwhite Ct, Suite 100, Boise, ID 83706
Telephone: (208) 342-5898 Ext. 115 Email: tsoppe@wegai.com
Supports and Attachments Preapproval
<ul> <li>Supports and attachments are preapproved under OPM- (Separate application for OSHPD Preapproval of Manufacturer's Certification (OPM) of Supports and attachments is required)</li> <li>Supports and attachments are not preapproved</li> </ul>
Certification Method
<ul> <li>Testing in accordance with: X ICC-ES AC156</li> <li>Other (Please Specify):</li></ul>
Testing Laboratory
Company Name: IABG mbH
Contact Name: Dr. Steffen Roedling
Mailing Address: Einsteinstrasse 20. Ottobrunn. Germany D-85521

 Telephone:
 +49 (0) 89 / 6088-2052
 Email:
 roedling@iabg.de

"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"



Page 2 of 3



Seismic Paramete	ers
------------------	-----

Design in accordance with ASCE 7-10 Chapter 13: $\square$ Yes $\square$ No Design Basis of Equipment or Components ( $F_p/W_p$ ) = Multiple, See attachment
$S_{DS}$ (Design spectral response acceleration at short period, g) = 2.2 for z/h = 1.0 and 2.5 for z/h = 0
$a_p$ (In-structure equipment or component amplification factor) = Multiple, See attachment
$R_{p}$ (Equipment or component response modification factor) = Multiple, See attachment
$\Omega_0$ (System overstrength factor) = <u>Multiple</u> , See attachment
$I_p$ (Importance factor) = 1.5
z/h (Height factor ratio) = <u>1.0 at S<sub>DS</sub> = 2.2g and 0 at S<sub>DS</sub> = 2.5g</u>
Equipment or Component Natural Frequencies (Hz) = <u>Multiple, See attachment</u>
Overall dimensions and weight (or range thereof) = <u>Multiple, See attachment</u>
Equipment or Components @ grade designed in accordance with ASCE 7-10 Chapter 15: 🗌 Yes 🛛 No
Design Basis of Equipment or Components (V/W) =
S <sub>DS</sub> (Design spectral response acceleration at short period, g) =
S <sub>D1</sub> (Design spectral response acceleration at 1 second period, g) =
R (Response modification coefficient ) =
$\Omega_0$ (System overstrength factor) =
Cd (Deflection amplification factor) =
$I_p$ (Importance factor) = 1.5
Height to Center of Gravity above base =
Equipment or Component Natural Frequencies (Hz) =
Overall dimensions and weight (or range thereof) =
Tank(s) designed in accordance with ASME BPVC, 2015: Yes No
List of Attachments Supporting Special Seismic Certification
☑ Test Report(s) □ Drawings □ Calculations ☑ Manufacturer's Catalog
Other(s) (Please Specify): Certified System Matrix, UUT Summary Sheets, Subcomponent Certification Letter
OSHPD Approval (For Office Use Only) – Approval Expires on December 31, 2022
Atte
Signature: Date: Date: May 17, 2017
Print Name: Ali Sumer Title: DSE
Special Seismic Certification Valid Up to : S <sub>DS</sub> (g) = <u>See Above</u> z/h = <u>See Above</u>
Condition of Approval (if applicable):
"Access to Safe, Quality Healthcare Environments that Meet California's Diverse and Dynamic Needs"
STATE OF CALIFORNIA – HEALTH AND HUMAN SERVICES AGENCY
OSH-FD-759 (REV 12/16/15) Page 3 of 3

#### SIEMENS HEALTHCARE GmbH SPECIAL SEISMIC CERTIFICATION CERTIFIED SYSTEM AND COMPONENTS



#### Manufacturer: Siemens Healthcare GmbH

#### System: Multitom Rax

System Component	Siemens	Ι	Dimensio	ns (in)	Weight	Mounting	UUT
System Component	Part Number	Width	Length	Height	(lb)	Mounting	001
Tube Ceiling Stand 3D V	07042075	44.8	34.3	33.5-118.6 <sup>2)</sup>	585	ceiling	UUT-1
Detector Ceiling Stand 3D V	07042026	32.4	72.6	41.3-124.2 <sup>2)</sup>	546	ceiling	UUT-2
RF Table / Tabletop	10092902 10882788	29.5	125.5	19.7-37.6 <sup>3)</sup>	668.8 <sup>3)</sup>	floor	UUT-3
Generator Polvdoros F80-2	10096925	31.5	17.1	86.7	862	floor/wall	UUT-4
PC (W550)	11105103	13.4	27.4	22.8	90	floor	UUT-6ab
FLC PC (W550 RAD)	11105102	13.4	27.4	22.8	87	floor	UUT-7ab

<sup>1)</sup> All components are manufactured by Siemens Healthcare GmbH unless noted. The part numbers listed uniquely identify the type of component, manufacturer, and material of construction for each sub-component within the tested units.

<sup>2)</sup> Tube and Detector were subjected to two tests: the first test was performed in the normal operating position of 63in and the second test was performed with the extended position of 90.5in (both measured to focal point).

<sup>3)</sup> Weight does not included simulated patient weight. The patient table was subjected to two tests; the first test was performed in the extended mid-possition 29.5in with 308lbs simulated patient weight and the second test was performed in the bottom possition 20.5in with 529lb simulated patient weight.

SEISMIC CERTIFICATION LIMITS										
System Component	Code	S <sub>DS</sub> (g)	z / h	I <sub>P</sub>	a <sub>P</sub>	R <sub>P</sub>	Ω <sub>0</sub>	$\mathbf{F}_{\mathbf{P}}$ / $\mathbf{W}_{\mathbf{P}}$		
Tube Ceiling Stand		2.2	1.0	1.50	2.5	2.5	2.0	3.96		
Tube Cennig Stand		2.5	0	1.50	2.3	2.3	2.0	1.50		
Detector Ceiling Stand		2.2	1.0	1.50	2.5	2.5	2.0	3.96		
Detector Cennig Stand		2.5	0	1.50	2.3	2.3	2.0	1.50		
RF Table / Tabletop	16 10	2.2	1.0	1.50	1.0 2.5	1.5 6.0	1.5	2.64		
KI Table / Tabletop	201 E7-1	2.5	0	1.50			1.3	1.13		
Generator Polydoros F80-2	CBC 201 ASCE7-1	2.2	1.0	1.50			2.0	1.65		
Generator Forvitoros 1780-2	A C	2.5	0	1.50				1.13		
PC (W550)		2.2	1.0	1.50	1.0	2.5	2.0	1.58		
rC (W330)		2.5	0	1.30	1.0	2.3	2.0	1.13		
FLC PC (W550 RAD)		2.2	1.0	1.50	1.0	2.5	2.0	1.58		
FLC FC (W 330  KAD)		2.5	0	1.30	1.0	2.5	2.0	1.13		

# UNIT UNDER TEST (UUT) SUMMARY SHEET



Mounting Details: Runner rails bolt to unistrut with 2 - M10 screws at each intersecting location

<image/>										
Manufacture	r: Siemens Health	care GmbH	I							
Component:	Tube Ceiling Stan	d					2075 / 100	1		
Component: UUT Functio	Tube Ceiling Stand	d used for ma	aking X-ray	/ exposur			12075 / 100	1		
Component: UUT Functio UUT Descrip	Tube Ceiling Stand n: Digital system tion: Component of	d used for ma of Multiton	aking X-ray n Rax Syste	/ exposur	es of the bo	ody		1		
Component: UUT Functio UUT Descrip	Tube Ceiling Stand	d used for ma of Multiton pratory, Ger	aking X-ray n Rax Syste rmany	/ exposur em	es of the bo			1		
Component: UUT Functio UUT Descrip	Tube Ceiling Stand n: Digital system tion: Component of	d used for ma of Multiton oratory, Ges	aking X-ray n Rax Systo rmany UUT PRO	/ exposur em <b>PERTIE</b>	es of the bo	ody : June 2016	5			
Component: UUT Functio UUT Descrip	Tube Ceiling Stand n: Digital system tion: Component on: IABG Test Labo	d used for ma of Multiton oratory, Gen Dimensio	aking X-ray n Rax Syste rmany <b>UUT PRO</b> ons (inches)	exposur em PERTIE	es of the bo Test Date S	ody : June 2016				
Component: UUT Functio UUT Descrip Test Location Weight (lb)	Tube Ceiling Stand n: Digital system ( tion: Component of n: IABG Test Labo Width <sup>1)</sup>	d used for ma of Multiton oratory, Gen Dimensio Dep	aking X-ray n Rax Syste rmany UUT PRO ons (inches) oth <sup>1)</sup>	em PERTIE Tested F	es of the bo	ody : June 2016 Natur	5 al Fequenc	y (Hz)		
Component: UUT Functio UUT Descrip Test Location	Tube Ceiling Stand n: Digital system tion: Component on: IABG Test Labo	d used for ma of Multiton oratory, Gen Dimensio Dep	aking X-ray n Rax Syste rmany <b>UUT PRO</b> ons (inches)	PERTIE Tested E	es of the bo Test Date S Extension <sup>2)</sup>	ody : June 2016 Natur FB	al Fequenc SS	y (Hz) V		
Component: UUT Functio UUT Descrip Test Location Weight (lb) 585	Tube Ceiling Stand n: Digital system ( tion: Component of n: IABG Test Labo Width <sup>1)</sup>	d used for ma of Multiton pratory, Gen Dimensio Dimensio Dep 34 tal directions contal placem the first test ad test was pe	aking X-ray n Rax Syste rmany UUT PRO ons (inches) oth <sup>1)</sup> .3" along the tra ent has no aff was performed with	PERTIE PERTIE Tested E 6 9 ck that mod fect on seis ed in the no the extend	es of the bo Test Date S Extension <sup>2)</sup> 3.0" 0.5" unts to the cei mic loading. ormal operation led position o	bdy June 2016 Natur FB -4.4 -3.4 iling adapter. ag extended p	5 al Fequency SS -3.7- -2.5- The system is osition of 63.0	y (Hz) V -4.9- -4.6- s uniformly Oin		
Component: UUT Functio UUT Descrip Test Location Weight (lb) 585 <sup>1)</sup> The UUT is op connected to the <sup>2)</sup> The UUT was (measured to foc	Tube Ceiling Stand n: Digital system of tion: Component of n: IABG Test Labo Width <sup>1)</sup> 44.8" berable in both horizon ceiling therefore horiz subjected to two tests: al point) and the secor	d used for ma of Multiton oratory, Gen Dimensio Dimensio Dep 34 tal directions contal placem the first test ad test was pe SEISM	aking X-ray n Rax Syste rmany UUT PRO ons (inches) oth <sup>1)</sup> .3" along the tra ent has no aff was performed with IIC TEST	PERTIE PERTIE Tested E 6 6 9 6 6 6 7 9 6 6 6 7 9 6 6 7 9 7 8 7 8 7 8 7 8 7 8 7 8 7 8 7 8 7 8	Test Date Test Date S Extension <sup>2)</sup> 3.0" 0.5" unts to the cei mic loading. ormal operation led position o ETERS	bdy June 2016 Natur FB -4.4 -3.4 iling adapter. ag extended p f 90.5in (mea	5 al Fequency SS -3.7- -2.5- The system is osition of 63.0 issured to focal	y (Hz) V -4.9- -4.6- s uniformly Oin point).		
Component: UUT Functio UUT Descrip Test Location Weight (lb) 585 <sup>1)</sup> The UUT is op connected to the <sup>2)</sup> The UUT was (measured to foc Building Co	Tube Ceiling Stand n: Digital system of tion: Component of n: IABG Test Labor Width <sup>1)</sup> 44.8" perable in both horizon ceiling therefore horiz subjected to two tests:	d used for ma of Multiton pratory, Gen Dimensio Dimensio Dep 34 tal directions contal placem the first test ad test was pe	aking X-ray n Rax Syste rmany UUT PRO ons (inches) oth <sup>1)</sup> .3" along the tra ent has no aff was performed with	PERTIE PERTIE Tested E 6 9 ck that mod fect on seis ed in the no the extend	es of the bo Test Date S Extension <sup>2)</sup> 3.0" 0.5" unts to the cei mic loading. ormal operation led position o	bdy June 2016 Natur FB -4.4 -3.4 iling adapter. ag extended p	5 al Fequency SS -3.7- -2.5- The system is osition of 63.0	y (Hz) V -4.9- -4.6- s uniformly Oin		

# UNIT UNDER TEST (UUT) SUMMARY SHEET



Mounting Details: Runner rails bolt to unistrut with 2 - M10 screws at each intersecting location

<image/>										
	er: Siemens Health Detector Ceiling S		1	Model /	Serial Nur	<b>nber:</b> 0704	2026 / 100	1		
_	<b>n:</b> Digital system		aking X-ray							
UUT Descrip	tion: Component	of Multiton	n Rax Syste	em		_				
	n: IABG Test Labo				Test Date:	June 2016	<u>5</u>			
		ו	UUT PRO	PERTIE	S					
W ' 1 ( /11 )		Dimensio	ns (inches)	)		Natural Fequency (Hz)				
Weight (lb)	Width <sup>1)</sup>	Dep	oth <sup>1)</sup>	Tested I	Extension <sup>2)</sup>	FB	SS	V		
546	32.4"		.6"	6	3.0"	4.3	-3.3-	<del>-3.1-</del>		
					0.5"	-3.4-	-2.6-	-3.7-		
<ul> <li><sup>1)</sup> The UUT is operable in both horizontal directions along the track that mounts to the ceiling adapter. The system is uniformly connected to the ceiling therefore horizontal placement has no affect on seismic loading.</li> <li><sup>2)</sup> The UUT was subjected to two tests: the first test was performed in the normal operating extended position of 63.0in (measured to focal point) and the second test was performed with the extended position of 90.5in (measured to focal point).</li> </ul>										
Duilding C.	da / Taat Oritaria			PAKAN I <sub>P</sub>		$A_{pro}$	<b>Α</b> <sub></sub> (σ)	$A_{pro}$		
Building Co	de / Test Criteria	$\frac{S_{DS}(g)}{2.20}$	z / h 1.0	п <sub>Р</sub> 1.5	A <sub>FLX-H</sub> (g) 3.52	A <sub>RIG-H</sub> (g) 2.64	$A_{FLX-V}(g)$	$A_{RIG-V}(g)$		
CBC 2016 /	ICC-ES AC156	2.20	0.0	1.5	5.52	2.04	1.67	0.67		

### UNIT UNDER TEST (UUT) SUMMARY SHEET



Mounting Details: Rigid Floor mounted with 6 - M12 bolts



Note: The unit was full of contents during testing and remained fuctional before and after the ICC-ES AC156 test. The unit maintained structural integrity during and after the ICC-ES AC156 Test.

# UNIT UNDER TEST (UUT) SUMMARY SHEET



Mounting Details: Floor/Wall mount with 4 - M12 bolts to floor and 2 - M6 bolts to wall fixture



#### UUT-6a

### UNIT UNDER TEST (UUT) SUMMARY SHEET



**Mounting Details:** Rigid Floor mounting using Siemens provided seismic restraint kit for Configuration A. Siesmic restraint kit includes two 1" wide hand tightened cam buckle straps (830lb WLL) looped thru angle brackets positioned on each side of the unit. The four angle brackets are attached to the table with individual M10 bolts.


<b>Component:</b>	PC (W550)			Model /	Serial Nun	nber: 1110	5103 / 1185	5	
<b>UUT Functio</b>	n: Computer for d	ata acquisit	ion, image	reconstr	uction, and	processing			
UUT Descrip	tion: Component of	of Multiton	n Rax Syste	em					
Test Location	n: IABG Test Labo	oratory, Gei	rmany		Test Date:	January 2	017		
		I	UUT PRO	PERTIE	S				
W/ ' 1 / /11 \		Dimensio	ns (inches)	)		Natural Fequenc		(Hz)	
Weight (lb)	Width	De			eight	FB	SS	V	
90	13.4"	27.	.4"	2	2.8"	24	18.6	>33	
		SEISM	IC TEST	PARAM	ETERS				
Building Co	de / Test Criteria	$S_{DS}(g)$	z / h	I <sub>P</sub>	$A_{FLX-H}(g)$	$A_{RIG-H}(g)$	$A_{FLX-V}\left(g ight)$	$A_{RIG-V}(g)$	
CBC 2016 /	ICC-ES AC156	2.20 2.50	1.0 0.0	1.5 1.5	3.52	2.64	1.67	0.67	
	as full of contents dur tural integrity during a			uctional be	fore and after	the ICC-ES	AC156 test. T	he unit	

#### UUT-6b

### UNIT UNDER TEST (UUT) SUMMARY SHEET



**Mounting Details:** Rigid Floor mounting using Siemens provided seismic restraint kit for Configuration B. Siesmic restraint kit includes one 1" wide hand tightened cam buckle strap (830lb WLL) looped thru angle brackets positioned on each side of the unit. Angle brackets are also positioned on the front and back of the unit. All four angle brackets are attached to the table with individual M10 bolts.



#### UUT-7a

### UNIT UNDER TEST (UUT) SUMMARY SHEET



**Mounting Details:** Rigid Floor mounting using Siemens provided seismic restraint kit for Configuration A. Siesmic restraint kit includes two 1" wide hand tightened cam buckle straps (830lb WLL) looped thru angle brackets positioned on each side of the unit. The four angle brackets are attached to the table with individual M10 bolts.

Manufacture	er: Siemens Health	care GmbH	I						
<b>Component:</b>	FLC PC (W550 R.	AD)		Model /	Serial Nur	<b>nber:</b> 1110	05102 / 116	2	
UUT Functio	<b>n:</b> Computer for d	ata acquisit	tion, image	reconstr	uction, and	processing			
<b>UUT Descrip</b>	tion: Component	of Multiton	n Rax Syste	em					
	n: IABG Test Labo		-		Test Date	January 2	017		
			UUT PRO	PERTIE	S	-			
			ns (inches)			Natur	latural Fequency (Hz)		
Weight (lb)	Width		pth		eight	FB	SS	V	
87	13.4"		.4"		2.8"	27.1	18.1	>33	
		SEISM	IC TEST	PARAM	IETERS				
Building Co	de / Test Criteria	$S_{DS}(g)$	z / h	I <sub>P</sub>	$A_{FLX-H}(g)$	$A_{RIG-H}(g)$	$A_{FLX-V}(g)$	$A_{RIG-V}(g)$	
	ICC-ES AC156	2.20 2.50	1.0	1.5 1.5	3.52	2.64	1.67	0.67	
	vas full of contents dur tural integrity during a	ing testing an	d remained f	uctional be	fore and after	the ICC-ES			

# UUT-7b

### UNIT UNDER TEST (UUT) SUMMARY SHEET



**Mounting Details:** Rigid Floor mounting using Siemens provided seismic restraint kit for Configuration B. Siesmic restraint kit includes one 1" wide hand tightened cam buckle strap (830lb WLL) looped thru angle brackets positioned on each side of the unit. Angle brackets are also positioned on the front and back of the unit. All four angle brackets are attached to the table with individual M10 bolts.

